Easy and simple measurement of salt % by conductivity method!

Up-to-the-minute funnel-shaped sensor section.

Compact body for two-way measurements being held by hand and placed on desk.

ES-421 SALT-METER

Sallmer 1971 50%

Scale of salt % : sodium chloride (g/100g)

..... Sensor section

The funnel-shaped sample receiver employed in the sensor section realizes measurement with a small amount of sample. Cleaning of the sensor section after measurement is very easy by just wiping it with tissues moistened with water. Differently from salt-meters employing glass electrodes, the ES-421 is free from fear of damage of the electrode and from special care to keep and store it.

Sample stage

The sample stage made of corrosion resisting hard plastics serves for a long time with cleanliness of the surface.

Little sample is enough for measurement!


Cleaning by just wipeover!

Simple calibration with airl


Easy and simple operation

The digital salt-meter ES-421 employs a funnel-shaped sensor section (conductivity method), which measures the salt % of a sample (dilute solution of soy sauce, ketchup, mayonnaise, Worcester sauce, squeezed juice of pickles, etc.) poured into it by a little amount (1.0ml approx. is enough for measurement).

*For measurement of a sample containing ingredients other than salt, dilute the sample so that its Brix index (percentage of soluble solid content) is less than 6%. For example, dilute soy sauce, Worcester sauce with distilled water 10 times (in weight) as thin as the original for measurement, and the actual salt % is obtained by multipliying the measured value by 10.



Apply drops of sample into the sensor section until it fills above the indication line.



Press the Start/Off switch.



After arrow flashes three times, the display indicates the salinity of the measured sample in percentage.

......... High reliable digital indicator

A large-sized LCD employed for indication of measured value enables the operator to check the test result easily and promptly with its excellent readability.

...... Excellent operability

Simple design with only two switches on the panel. Turning power on/off, zero setting, start of measurement can be operated with one touch of a switch.

Moreover, two switches designed considerably big for easy pressing facilitate correct and quick measurement for every operator. The ES-421 indicates salt % of sample in digits just by filling up the funnel-shaped sensor section with a small amount of sample (1.0ml approx.) and pressing the Start/Off switch.

•Note:

For measurement of food containing high conductive ingredients other than salt, it may be required to dilute sample specially or to prepare a special conversion table.